## **UC-03 Sequence-Diagram-XDS**

Scenario: Patient Portal Retrieves PS from XDS Document Registry – (XDS IHE Profile)

Assumption: Patient Summary is stored in Local (Decentralized) Document Repository (XDS style repository)

Alternative Option: This sequence diagram provides the option of using the XDS IHE profile to give jurisdictions and vendors option to leverage existing XDS-based implementations in the ecosystem.

Note: This sequence diagram includes a custom Canadian IHE profile (e.g. CA:FMT) that handles transformations to and from various formats (e.g. FHIR to PDF, CDA, etc.). Additional details will be included in the PS-CA Interoperability Specifications.


## Sequence Diagram Overview:

Below provides guidance on how to read the sequence diagram:

- This sequence diagram illustrates how the different standardized actors of system should interact with each other to carry out specific standardized transactions, and the order in which the transactions and interactions occur when Use Case 1 of the Patient Summary-CA is executed.
- The legend on the bottom right corner describes the different system components, actors and transactions that are necessary to carry out this particular use case.
- The green swim lane is a simplified view of the actors and transactions required by the Foundational IHE Profiles, in addition to the other ones
  that are not explicitly illustrated on the diagram (e.g. ATNA, CT) but included as a white note. These are pre-requisite conditions for this particular
  use case and it is assumed that these will be satisfied.
- The blue swim lanes groups sequence of processes (along with their required actors and transactions) that are needed to occur to satisfy this
  particular use case. These are to be read from left to right and top to bottom.
- The red note boxes describe important call outs, information and notes that provide more context for the sequence diagram.
- This sequence diagram includes a custom Canadian IHE profile (e.g. CA:FMT) that handles transformations to and from various formats (e.g.
  FHIR to PDF, CDA, etc.). Additional details will be included in the PS-CA Interoperability Specifications.
- More information about those details of the Foundational IHE Profiles an be found here.

## UC-03 Patient Views/Obtains Personal PS-CA - XDS

Patient Portal Retrieves PS from XDS Document Registry – (XDS IHE Profile)
PS is stored in Local (Decentralized) Document Repository (XDS style repository)





Client is logged into the system (IUA\*)
Client obtained a valid access token
from the Authorization Server that is
used with each transaction (IUA\*)
All communication is done through
secure channels (ATNA)
System time is synchronized among all
components (CT)

